


Product Range

Environmental Monitoring Solutions



Simple Color Tests

(Manual Input)




Residue Test

- SpotCheck Plus**
Glucose & Lactose
- PRO-Clean**
Protein
- AllerSnap**
High-Sensitivity, Broad Spectrum, Protein
- Allerflow | Gluten**
High-Sensitivity, Gluten Specific


Pathogen Screening Test

- InSite® | Listeria**
Ready-to-Eat, Dairy, Produce, Meat
- InSite® | Salmonella**
Ready-to-Eat, Dairy, Produce, Meat
- InSite® | L. mono Glo**
Easy to interpret color change



Secure data integration
from multiple test type
and data sources

EnSURE Touch
with Charging Dock



CalCheck



Calibrates your luminometer in less than 1 minute



Luminometer Tests

ATP Bioluminescence
(Automatic Syncing Input)

ATP Cleaning Verification Test

- UltraSnap**
Surface ATP Test
- SuperSnap**
High-Sensitivity Surface ATP Test
- AquaSnap TOTAL**
Water ATP Test
- AquaSnap FREE**
Water ATP Test

Enzyme and Efficiency Test

- ZymoSnap ALP**
Alkaline Phosphatase, Dairy
- CROSSCHECK**
Acid Phosphatase, Raw Meat

Indicator Organism Test

MicroSnap®

MicroSnap TOTAL

Enumerate - 7 hrs

MicroSnap EB

Enumerate - 6 hrs | Presence/Absence - 8 hrs

MicroSnap COLIFORM

Enumerate - 6 hrs | Presence/Absence - 8 hrs

MicroSnap E. coli

Enumerate - 6 - 7 hrs | Presence/Absence - 8 hrs